

Peak Analysis Tool User Manual

Contents

Peak Analysis Tool User Manual	1
Table of Contents	1
Introduction	1
Getting Started	2
Workflow Overview	2
Loading Data	2
Process & Detect	2
Analysis & Export	3
Double Peak Analysis	3
Exporting Results	3
Preferences	3

Peak Analysis Tool User Manual

Table of Contents

1. Introduction
2. Getting Started
3. Workflow Overview
4. Loading Data
5. Processing & Detection
6. Analysis & Visualization
7. Double Peak Analysis
8. Exporting Results
9. Preferences

Introduction

The Peak Analysis Tool is a high-precision application for detecting, analyzing, and characterizing peaks in time-series data. It is designed for scientific use cases, such as single-molecule fluorescence or particle counting, where accurate peak detection and noise filtering are critical.

Getting Started

Installation

- **Windows:** Run `PeakService.exe`.
- **Source:** Run `python -m service.app`.

The application interface will open automatically in your default web browser at `http://localhost:8765`.

Workflow Overview

The application follows a structured workflow: 1. **Load:** Import raw data files. 2. **Process & Detect:** Filter noise and identify peaks. 3. **Analyze:** View time-series results and distributions. 4. **Export:** Save data and plots.

Loading Data

Navigate to the **Load Data** page to begin.

- **File Selection:** Click “Select Files” to choose `.txt`, `.xls`, or `.xlsx` files.
- **Protocol Information:** Enter metadata about your experiment (Date, Setup, Sample ID, Laser Power, etc.) to keep it associated with the dataset.
- **Settings:**
 - **Time Resolution:** Set the time per sample (e.g., 0.1 ms).
 - **Photon Correction:** Enable dead-time correction if using a photon counter. Enter the dead time in nanoseconds.

Recent sessions are saved and can be reloaded with a single click from the “Recent Sessions” list.

Process & Detect

This is the core analysis step.

1. Signal Filtering

Apply a digital filter to reduce noise before detection. - **Filter Type:** - **Butterworth:** General-purpose noise reduction. - **Savitzky-Golay:** Preserves peak shape better while smoothing. - **Parameters:** - **Cutoff Frequency (Hz):** Frequencies above this are removed. - **Order:** Sharpness of the filter cutoff. - **Window Length:** Number of points for smoothing (Savitzky-Golay).

2. Peak Detection

Configure algorithms to find valid peaks. - **Prominence:** How much a peak stands out from the surrounding baseline. - **Min Distance:** Minimum samples between peaks. - **Width Min/Max:** Accepted width range in milliseconds.

- **Rel Height:** The relative height (0.0-1.0) at which width is measured.
- **Prominence Ratio:** Filters noise by comparing prominence to absolute amplitude.

Click **Detect Peaks** to run the algorithm. Use **Auto Threshold** for an initial estimate.

Analysis & Export

View the detected peaks in context of the full time series.

- **Interactive Plot:** Zoom and pan to inspect individual peaks.
- **Throughput:** See the rate of peaks over time.
- **Histograms:** Analyze distributions of peak amplitudes and intervals.
- **Export:** Use the integrated Export controls to save peak data (CSV) or download the current chart view as an image.

Double Peak Analysis

For experiments involving paired events (e.g., double-labeled particles).

- **Pairing Logic:** Peaks are paired based on their temporal proximity.
- **Metrics:**
 - **Distance:** Time between paired peaks.
 - **Prominence Ratio:** Ratio of the second peak's intensity to the first.
 - **Width Ratio:** Comparison of peak shapes.
- **Thresholds:** Set min/max ranges for these metrics to filter valid pairs.

Exporting Results

Save your work for publication or further analysis.

- **Data Formats:** CSV, Excel, Text.
- **Metadata:** Option to include protocol details in the header.
- **Images:** Click the **Camera** icon on any chart to save a high-resolution PNG.

Preferences

Customize the application appearance (Light/Dark mode) and default analysis parameters.